

Intended Use

The purpose of the **MBRAUN** gas purification system, together with a glove box, is for the enabling and maintaining a pure inert gas atmosphere inside a hermetically sealed enclosure. Materials that are sensitive to moisture and/or oxygen are handled by using the attached gloves or additionally specially designed handling systems. The system is intended for professional use only.

NOTICE

The system will require certain modification for working with delicate or dangerous materials.

Examples of such cases include:

- pharmaceutical or nuclear applications
- working with substances that will lead to dangerous situations if exposed to air in case of a failure of the system
- very expensive materials that might be destroyed if exposed to air in case of a failure of the system

These application types must be discussed with **MBRAUN** prior to operating the system.