

— *turning ideas into solutions*

**MBRAUN**

## **MB GS-35** REFRIGERATOR



- Variable shelves
- External compressor
- Stainless Steel 1.4301 (US type 304)
- PLC temperature control

**INERTGAS TECHNOLOGY**



## MB GS-35 TECHNICAL DATA

Temperature *range*	From +10°C to -35°C
Hysteresis	+/- 3°K
Inside dimensions	423 mm (h) x 266 mm (w) x 162.5 mm (d)
Volume	18 liters
Body and door	brushed Stainless Steel, type 1.4301
Shelves	1 fixed, 2 variable shelves on 5 possible levels, Stainless Steel, type 1.4301 (US type 304). The shelves are secured against dumps and falling out.
Coolant	R404A
Electrical power	100 - 230 V, 50 - 60 Hz, 5 A

\* A retrofit version with Beta controller has temperature range of -10°C to -35°C and hysteresis +/- 5°K is available

**M. Braun Inertgas-Systeme GmbH (Headquarters)**  
Dieselstr. 31 · 85748 Garching · Germany  
Phone: +49 89 32669-0 · Fax: +49 89 32669-105  
Web: [www.mbraun.com](http://www.mbraun.com) · E-mail: [info@mbraun.de](mailto:info@mbraun.de)

**M. Braun Incorporated**  
14 Marin Way · Stratham, NH · 03885 · USA  
Phone: +1 603 773-9333 · Fax: +1 603 773-0008  
Web: [www.mbraun.com](http://www.mbraun.com) · E-mail: [info@mbraunusa.com](mailto:info@mbraunusa.com)

MB GS 35-EZ012005 · © M. Braun GmbH · Subject to be changed without notice!