

MB-GS-35

FREEZER



- Adjustable temperature range from -35 to +10°C
- Variable shelf levels
- Energy efficient external compressor cooling
- Side or rear wall mounted version
- Resistant to most commonly used solvents
- Ergonomic design allows accessing the entire internal space

MB-GS-35

FREEZER



Clean. Engineering. Expertise.



TECHNICAL DATA	
Position	Left, right or rear wall
Temperature range	From +10°C to -35°C
Hysteresis	3°K
Inside dimensions	423 mm (h) x 262 mm (w) x 162 mm (d)
Volume	18 l
Body and door	Brushed Stainless Steel, type 1.4301
Shelves	3 variable shelves on 5 possible levels, stainless steel, type 1.4301 (US type 304). The shelves are secured against dumps and falling out
Ambient temperature range	0°C to +35°C
Coolant	R404a
Electrical power	115 - 230 V, 50 - 60 Hz, 5 A

Technical note: Dependent upon operating conditions.

M. Braun Inertgas-Systeme GmbH (Headquarters)
Dieselstr. 31 • D-85748 Garching • Germany
Phone: +49 89 32669-0 • Fax: +49 89 32669-105
Web: www.mbraun.de
E-Mail Sales: info@mbraun.de
E-Mail Service: service@mbraun.de

M. Braun Incorporated
14 Marin Way • Stratham, NH • 03885 • USA
Phone: +1 (603) 773 9333 • Fax: +1 (603) 773 0008
Web: www.mbraunusa.com
E-Mail Sales: info@mbraunusa.com
E-Mail Service: service@mbraunusa.com